

life - P  INT



**BIPHASIC
DEFIBRILLATOR**


0535

w w w . m e t s i s m e d i c a l . c o m

PRO BIPHASIC 270 JOULE

Technical Specifications

Wave Form	Biphasic
Operating Mode	Synchronous / Asynchronous
Energy Selection	2 - 270 Joule
ECG Derivations	DI-DII-DIII-AvR-AvL-AvF-C1...C6
ECG Filters	Mon1-Mon2-Strong-Surgery-Diagnostic
Heart Rate Alarm Limits	15 - 300
Freeze	Yes
Monitor Type	LCD
Monitor Resolution	320x240
Monitor Size	5.5 inch
Monitor Speed	25mm/sec
Printer Type	Thermal Printer
Printer Speed	25-50mm/sec
Paper Size	50mm x 30m
Max. Discharging	100+
ECG Monitoring Time	4 Hours
Working Conditions	0/40 °C ve %30-95Rh
Storing Conditions	-20 / +55 °C
Dimensions	403x152x324 mm
Weight	8,5 Kg
Class	II B
Main Supply Voltage	AC230V/1A 50-60 Hz
Direct Working From Main Voltage	Optional
Pulseoximeter	Optional
Standart Accessories	ECG Cable and Paper ECG Electrode - Power Cable

Opsiyonel Aksesuarlar



Stand



Trolley



ANKARA CYBERPARK
To serve you better, we continue to work in one of the best "technology development zone" of Turkey



MANUFACTURER

Head Office: İvedik Organize Sanayi Bölgesi Sogutmacilar Sitesi 655. Sok.No:11 06370 Yenimahalle/ANKARA
Tel: +90 (312) 3957770 **Fax:** +90 (312) 3957772 **E-Mail:** info@metsismedikal.com
Branch Office: Cyberpark Ankara Teknoloji Gelistirme Bölgesi Cyberplaza B-Blok 2.Kat No:B210 Bilkent/ANKARA

