




		FETAL DOPPLERS		
	Sonotrax -A	<ul style="list-style-type: none">• sonotrax handheld Doppler can meet routine examination requirement of obstetricians, midwives, and pregnant women. It has manual heart rate processing mode. It has audio output, and can connect with earphone or recorder with audio input. Sonotrax has a slim, power saving and interchangeable ultrasound probe.• No display		
	Soline C	<ul style="list-style-type: none">• 45*825mm LCD display• FHR Performance• FHR Range: 50-240BPM.• Probe Frequency: 2.0MHz• AAA 1.5 V batteries• Working Frequency: 2.0MHz 10%• Resolution: 1BPM• Accuracy: 2BPM• Power consumption: <0.8W• Auto Shut-OFF: After 1 minute no signal, power off automatically.		