

GD400

Electric surgical unit

Guoteng™

**One Of The Professional Chinese
MANUFACTURERS
Special In
Medical Equipment**



All automatic digital display system of output watt facilities confirmation of user power.
Equipped with automatic visual & acoustic monitoring & warning system,
the power output is automatically cut off in any Interruption
Cutting power and coagulating power do not work at the same time.
This can be hand or foot-operated.

High frequency power is generated at the active electrode only,
to obviate inadvertent diathermy.
Seven working modes.

Double frequency output, wide range and constant power.

Blend modes are designated low and high.

300W power output, totally adapt to surgical operation.



Guoteng Science & Technology Development Co., Ltd. Zhuhai

Address: 6/F., Three-Dimension Science & Technology Building, No. 184 Baillan Road, Jida, Zhuhai, Guangdong, 519015 China
Tel: +86-756-3321860 Fax: +86-756-6181801 E-mail: market@guoteng.gd.cn

<http://www.guoteng.gd.cn>