

life - POINT



**WITHOUT MONITOR
DEFIBRILLATOR**

CE
0535

www.metsismedical.com

WITHOUT MONITOR 360 JOULE

Technical Specifications

Wave Form
Operating Mode
Energy Selection
Max. Discharging
Battery Charging Time
Working Conditions
Storing Conditions
Weight
Dimensions
Class
Main Supply Voltage
Direct Working From Main Voltage
Standart Accessories

Monophasic
Synchronous / Asynchronous
5 - 360 Joule
100+
5 Hours
0/40 °C and %30-95Rh
-20 / +55 °C
8 Kg
403x152x324 mm
II B
AC230V/1A 50-60 Hz
Optional
Power Cable

Optional Accessories



Stand



Trolley



ANKARA CYBERPARK
 To serve you better, we continue to work in one of the best "technology development zone" of Turkey



MANUFACTURER

Head Office: Ivedik Organize Sanayi Bolgesi Sogutmacilar Sitesi 655. Sok. No:11 06370 Yenimahalle/ANKARA
Tel: +90 (312) 3957770 **Fax:** +90 (312) 3957772 **E-Mail:** info@metsismedikal.com
Branch Office: Cyberpark Ankara Teknoloji Gelistirme Bolgesi Cyberplaza B-Blok 2.Kat No:B210 Bilkent/ANKARA

