

# AED 501 Life Mode



Model	AED 501
Product Type	Low-output CPR Machine (Folded)
ECG Analysis	8.5 Sec.
Hi-Cap Charging	8 Sec.
Output	50J~150J offered
Data Storage	100 (ECG storage)
Data Transfer Protocol	Bluetooth 2.0
Control Button	Power, Shock, Adults/kids
Battery Type	Li-Mn (12V DC 4.2A)
Shock Frequency	200, at least
Pad Type	Electrode Pad, disposable and adhesive
Dimension	230(W) x 320(L) x 85(H)mm
Weight	1.9kg
Warranty	5 Yrs. (2 Yrs. for Pads/4 Yrs. for Batteries)