

H6 Patient Monitor

Features

- > 8.4" Colour TFT Display, led backlight , 800x600 pixels.
- > Operating Mode Selection: Neonatal/Pediatric/Adult.
- > Visual and Audio alarms, adjustable.
- > Patient management, name and age and id.
- > Network with Central Station Software,
- > Built-in removable and rechargeable lithium battery,
- > easy operating UI with both Keypad menu.
- > Powerful data storage and management, Trend graph, Trend table, etc
- > Arrhythmia analysis and S-T segment analysis.
- > Compatible to Neonatal, Pediatric, Adult Patients

Standard parameters:

3/5-Lead ECG, Hwatime SpO2 , NIBP, RESP, 2-Temp, PR

Optional:

EtCO2 , Touchscreen , Thermal Recorder ,WLAN accessory , Nellcor-SPO2 ,2-IBP .

H6 Patient Monitor



H6 Patient Monitor

Hwatime

Hwatime Medical • Make your life better

7/F, Block A, Bofook Jewelry Industry Park, No.33,
Bulan Rd, Longgang District, Shenzhen City
Tel: 0086-755-28260108 Fax: 0086-755-85221277
E-mail: sales@hwatime.com



Handset two-dimensional code

website: www.hwatime.com

Customer service hotline : 4000 670 888

· Hwatime. All Rights Reserved. The product specifications are subject to change without further notice.

· See instructions for contraindication or precautions

H6 Patient Monitor

Economic, Durable design

ECG

Lead Mode 3 / 5 leads
LEAD Selection I,II, III, avR, avF, V
GAIN Manual 2.5mm/mv, 5, 10 & 20mm/mv
Auto HR & alarm Adult 15-300BPM, Nep/Ped 15-350BPM
Accuracy +/- 1% or +/- 1bpm
Resolution 1 BPM
Sweep speed 12.5mm, 25mm, 50mm /second
Frequency Response Surgery : 1 – 20 Hz,
Monitor : 0.5 – 40 Hz,
Diagnostic: 0.05 – 130 Hz
Arrhythmia analysis 13 Kinds

NIBP

Method Oscillometric
Mode Manual, Auto, Start
Interval in auto mode
1,2,3,4,5,10,15,30,60,90,120,180,240 & 480mints
Accuracy Max Mean Error : +/- 5mmHg
Alarm type
SYS, DIAS & MEAN, : SYS: 40-275 mmHg,
DIAS:10-215 mmHg, Mean :20- 165 mmHg

SpO2

Measuring Range 1-100%
ALARM RANGE 0-100%
Resolution 1%
AVERAGING 4, 8 & 16 seconds
ALARM DISPLAY 10 seconds
PULSE RATE & ALARM 25-240BPM
RESOLUTION 1 BPM
ACCURACY +/- 2BPM
+/- 2% - 70% - 100%

Pulse Rate

Source User Slectable: Smart, ECG, Pleth, NIBP
Range Adult 15-300 bpm
Pediatric 15-350 bpm
Accuracy ± 1 bpm or $\pm 1\%$ whichever is greater

Temperature

Channels 2
Range, Accuracy 0°-50°C : ± 0.1 °C
Display Resolution ± 0.1 °C
Proble YSI 400 and YSI 700 Series

Respiration Rate

Rate Adult: 0 ~ 120 rpm
Neo/Ped 0 ~ 150 rpm
Resolution 1rpm
Accuracy ± 2 breath/min

Printer (optional)

Type Thermal printer
Paper Speed 25 or 50 mm/sec

Trends

Memory
96 hr at 30 second intervals for a NIBP reading every 5 minutes
Display tabular; graphical

Physically

Power Source 110V-240V, 50/60Hz
Weight Approximate 4.2Kg
Audio volume adjustable 1 ~ 4 level
SCREEN SELECTION
standard, Oxy CRG, Big Font & Trend



SpO2 sensor



ECG cable



NIBP tube and cuff



Durable probes connector
with high quality



TEMP probe