

RT2-04-07 Vacuum Suction with Extraction cups

High Vacuum Extractor

SUCTION YX932D + silicone cups

(DISPOSABLE LINERS)

- South African proofed reliability
- Max.vacuum: 0.09MPa (680mmHg)
- Adjustable vacuum Range: 0.02MPa - Max.Vacuum
- Flow rate: 32L/Min
- Storage bottle: 2000ml x2
- Silicone cups sizes 1 x 40mm , 2 x 50mm ,
2 x 60mm
- Foot activation pedal

