

**Hwatime 华腾**

Hwatime Medical. Make life more brilliant



## i12/i10 Patient Monitor



The model of basic health care



## **i12/i10** Patient Monitor

**Hwatime**  
Hwatime Medical. Make life more brilliant

7/F, Block A, Bofook Jewelry Industry Park, No.33,  
Bulan Rd, Longgang District, Shenzhen City  
Tel: 0086-755-28260108 Fax: 0086-755-85221277  
E-mail: sales@hwatime.com



Handset two-dimensional code

Please log in for more details [www.hwatime.com](http://www.hwatime.com)

Customer service hotline : 400-8972-888

- Hwatime. All Rights Reserved. The product specifications are subject to change without further notice.
- See instructions for contraindication or precautions

# i12/i10 Patient Monitor



## Smart, Light, Easy operating

Suitable for bedside, ICU, OR



### Product feature

- Standard parameters: ECG, NIBP, RESP, PR, SpO2, TEMP
- Colorful and Clear 12.1" & 10.4" TFT screen, Backlight buttons
- Multiple display modes optional: Standard interface, Large font, ECG standard full display, OXY, Trend table, BP trend, Viewbed
- Ambulatory blood pressure technology, anti-movement
- Special design against high frequency surgical unit, and defibrillation protection
- Support Masimo SpO2-authorization partner
- 13 kinds of arrhythmia analysis
- 15 kinds of drug dose calculation
- Various languages operating systems
- Built-in detachable rechargeable lithium battery 4 hours battery life
- To connect CMS, other bed observation and software updating with wireless and wired mode

### Data & Storage – stable and fast

- 8000 groups NIBP measurements
- 680 hours Trend data and Trend graphs
- 200 groups Alarm events reviewing
- 2 hours Wave forms reviewing

### Alarming – safe and reliable

- 3 level audible and visual alarming
- Dual alarm light for physiological and technical alarming



### Standard interface diagram

Clear, st analysis, providing more comprehensive diagnosis information



### Full lead interface

7 lead/12 lead real time same-screen display, be convenient for medical staff to observe patient's dynamics in whole period



### Tendency diagram

According to the requirements of clinics, capture data of different period for retrospective analysis



### Big font interface

Be convenient for medical staff to observe patient's conditions from a long distance

Optional: EtCO2, Dual-IBP, 12-leads ECG, Touch screen, Built-in thermal printer



**Double invasive blood pressure (IBP)**  
Adopting double-channel measurement, and realizing continuous testing blood pressure variation of severe patients



**COM port**  
Support HD, VGA, USB, Internet connector output

### i12/i10 Patient Monitor

Suitable for Adult, Pediatric and Neonatal. Reliable quality gives confidence, stable system makes perfect. Powerful storage capacity ensures data safety.

Hwatile Medical regards quality as the core, relies on creative R&D, produces human friendly products, and doctors' reliable partner.

### Movable handle design

Lightweight portable, practical all-around, convenient operation

