



BN-100(Luxurious) Infant radiant warmer

Function features

Microprocessor based servo controlled temperature system Three modes: pre-warm, manual or automatic

Set temperature and skin temperature display separately

Special design inclined heater head is more convenience for X rays Carbon steel heater for more effectivity

Audible and visual Alarm functions :

Power failure, temperature deviation, temperature sensor failure, skin temperature failure, over temperature APGAR timer to record therapy time

Convenient X-ray cassette under infant bed Electric inclination of infant bed

LED observation light adjustable,brightness adjustable

Four Transparent protector folded with 2 infusion seals

LED phototherapy unit,Oxygen supply system

Bigger infant bed

Built in low pressure aspirator

With RS-232 connector,2 trays and 2 drawers

Optional functions



LCD Touch screen



Air temp&count down timer



Electronic Elevator



Air&Oxygen Blender



Gel mattress



Infant Resuscitator