



BN-100(Top Grade) Infant radiant warmer

Function features

Microprocessor based servo controlled temperature system

The microprocessor-controlled unit allows for three modes: pre-warm, manual and automatic

APGAR timer to record therapy time

Audible and visual Alarm functions :

Power failure, temperature deviation, temperature sensor failure, skin temperature failure, over temperature

LED observation light adjustable, brightness adjustable

Heater head can be adjusted $\pm 90^\circ$ horizontally

Convenient X-ray cassette under infant bed

Tilt angle of infant bed is adjustable

Four Transparent protector folded with 2 infusion seals

Integrated low pressure aspirator

LED phototherapy unit

PVC resuscitator

With RS-232 connector, 2 trays and 2 drawers