



BN-100(Standard) Infant radiant warmer

Function features

Microprocessor based servo controlled temperature system

The microprocessor-controlled unit allows for three modes :pre-warm ,manual or automatic

Timer to record therapy time

Audible and visual Alarm functions :power failure, temperature deviation,
temperature sensor failure , skin temperature failure,over temperature.

LED observation light adjustable,brightness adjustable

Heater head can be adjusted $\pm 90^\circ$ horizontally

Convenient X-ray cassette board under infant bed

Tilt angle of infant bed is adjustable

Four Transparent protector folded with 2 infusion seals,

With RS-232 connector ,2 tray and 2 drawers




LCD Touch screen



Big infant bed



Phototherapy BL-60D



Air&Oxygen Blender



Low Pressure Aspirator



Infant Resuscitator



Oxygen supply system



Electronic Elevator



Air temp&count down timer

Optional functions