



BN-100A(Standard)

Infant radiant warmer

Function features

Microprocessor based servo controlled temperature system

The microprocessor-controlled unit allows for two modes: manual and automatic

Set temperature and skin temperature are displayed by LED

Alarm functions :

Power failure, temperature deviation, temperature sensor failure, skin temperature failure, over temperature A second thermal cut-out function for more safety

Adopts anti-explosion micro crystal quartz infrared radiation tube

Four Transparent protector folded with 2 infusion seals

RS-232 connector, 1tray




Heating head with metal grid



Siderail colume



Aluminium alloy material



2 infusion seals