



BN-200Plus

Microprocessor based servo controlled temperature system

Three modes: pre-warm, manual and automatic

Audible and visual Alarm functions :

Power failure, temperature deviation, temperature sensor failure,
skin temperature failure, over temperature

Skin temperature, auxiliary temperature and ambient temperature are displayed separately

Built-in breath combination mode (negative pressure suction, continuous oxygen infusion, cardiopulmonary resuscitation)

10"LCD colorful touch screen

Carbon steel heater more effective

Separate keyboard lock to prevent misoperation

Special design inclined heater head is more convenience for X-ray image

Convenient X-ray cassette tray under infant bed

Integral timer to record therapy time

Electric infant bed step-less adjustable tilt $\pm 12^\circ$

Four transparent protector with 12 infusion seals can be
flipped and removed, damping slow down

The height of whole unit can be adjusted, pedal switch

LED observation light, lighting direction 360adjustable

The tray is adjustable up and down

With oxygen supply system, LED phototherapy,

baby weight system, pulse oximeter

With RS-232 connector and 1 big drawer